

RESPIRATORY DEVICES

PRODUCT	COMPANY	VALVE	MASK*	FLOW SIGNAL	CHAMBER VOLUME	NOTES	MORE INFORMATION
HAND-HELD SPACERS AND HOLDING CHAMBERS							
SPACERS							
OptiHaler	<i>Respironics</i>	none	none	none	60mL	<ul style="list-style-type: none"> • May be used as a carrying case for one MDI • Mouthpiece cap converts to finger grip to aid actuation • Latex-free 	(800) 345-6443 www.healthcare.philips.com/us_en/homehealth/respironics.wpd
PocketSpacer	<i>nSpire Health</i>	none	none	audio	110mL	<ul style="list-style-type: none"> • Pocket size • Latex-free 	(800) 205-7187 www.nspirehealth.com
			infant				
			small				
			medium				
			large				
VALVED HOLDING CHAMBER							
AeroChamber Plus Flow-Vu	<i>Forest</i>	one-way	none	audio + visual	149mL	<ul style="list-style-type: none"> • Delivers more medication than original AeroChamber • Anti-static valved holding chamber • Low effort valves • Latex-free 	(800) 678-1605 www.frx.com/products/AeroChamber.aspx
			small	visual			
			medium				
			large	audio + visual			
EasiVent	<i>Dey</i>	one-way	none	audio	137mL	<ul style="list-style-type: none"> • Latex-free 	(800) 755-5560 www.mylanspecialty.com/products/
PocketChamber	<i>nSpire Health</i>	one-way	none	audio	108mL	<ul style="list-style-type: none"> • Pocket size • Latex-free 	(800) 205-7187 www.nspirehealth.com
			infant				
			small				
			medium				
			large				
OptiChamber	<i>Respironics</i>	one-way	none	audio	218mL	<ul style="list-style-type: none"> • Maybe used as carrying case for one MDI • Removeable valve • Latex-free 	(800) 345-6443 www.healthcare.philips.com/us_en/homehealth/respironics.wpd
Vortex	<i>Pari</i>	one-way	none	visual	194mL	<ul style="list-style-type: none"> • Non-electrostatic metal chamber improves dose delivery consistency • Latex-free • Dishwasher safe • Compact • Universal receptacle 	(800) 327-8632 www.pari.com
			infant				
			pediatric				
			adult, regular				
			adult, large				
RESERVOIR BAG							
E-Z Spacer	<i>FSC Laboratories</i>	none	none	visual	700mL	<ul style="list-style-type: none"> • Do not use abrasive cleansers or dishwasher • Do not use if bag is torn • Collapses to small size for carrying • Bag is reusable 	(914) 251-5160 www.e-zspacer.com
InspirEase	<i>Schering</i>	none	none	audio	700mL	<ul style="list-style-type: none"> • Do not use abrasive cleansers or dishwasher • Do not use if bag is torn • Includes actuation aid • Replace bag every 2 to 3 weeks 	(800) 526-4099

NOTES

*none = available without mask. Suggested ages given; all masks should be measured according to manufacturer's recommendations and fit to form a tight seal with face. (Rev. 12/2012)